



The Chemical Company

Content  
Uniformity

**Lytrol® micro**

The new measure for bioavailability.




## The Preface

*In recent years an increasing number of poorly soluble substances has been developed.*



This requires the addition of solubilizers in order to make those substances bioavailable. As solid oral dosage forms are the preferred formulation, solubilizers are required with special galenical properties.

Poloxamers are compounds with a long history of use in the field of solubilization that have excellent properties. We are now able to offer micronized grades of these monographed compounds, which results in several advantages:




*Trademarks are owned  
by BASF Aktiengesellschaft*

- Effective solubilization and increased bioavailability, elimination of dose dumping
- Low toxicity
- Perfect content uniformity
- Average particle size within the same range as for active ingredients, which results in homogenous mixtures for granulation or tableting

- Robust production process  
Compatibility of the active ingredient and excipients
- No need for micronization of active ingredients for better solubility







***Lutrol® F***  
**+ active**



***Lutrol® micro***  
**+ active**

## The Product

For **Lutrol®micro 68**:

*a* = approx. 79

and *b* = approx. 28

The proportion of polyoxyethylene by weight is approx. 70 %


For **Lutrol®micro 127**:

*a* = approx. 101

and *b* = approx. 56

The proportion of polyoxyethylene by weight is approx. 80 %

The **Lutrol®micro** grades are micronized poloxamers working as solubilizing agents. A special micronization process is applied to achieve superior properties.



**Structural formula:** **Lutrol®micro 68** and **127** are polyoxyethylene-polyoxypropylene block copolymers with an average molecular weight of around 8,600 Daltons and 12,200 Daltons, respectively.

### Product form:

At room temperature **Lutrol®micro 68** and **127** appear as white powder prepared by micronization of prilled poloxamer grades. The products are manufactured by alkali-catalyzed polymerization of propylene oxide followed by ethoxylation of the polyoxypropylene unit with subsequent neutralization. The products contain butylhydroxytoluene as an antioxidant and have a weak odor.

### Pharmacopoeia:

Both products meet the requirements of the current USP-NF "Poloxamer", JPE "Polyoxyethylene/Polyoxypropylene Glycol" and Ph. Eur. family monograph "Poloxamers".



## Chemical nature:

The polymers are block copolymers of polyoxyethylene-polyoxypropylene. Both products have been micronized to an average particle size of approx. 50 µm.

	Lutrol®micro 68	Lutrol®micro 127
Molecular weight	7680–9510 g/mol	9840–14600 g/mol
pH-value (2.5 % in water)	5.0–7.5	5.0–7.5
Unsaturation	0.018–0.034 meq/g	0.031–0.065 meq/g
Melting point	52–57 °C	50–56 °C
Polyoxyethylene by weight	79.9–83.7 %	71.5–74.9 %
Heavy metals	max 10 ppm	max 10 ppm
Ethylene oxide	max 1 ppm	max 1 ppm
Propylene oxide	max 5 ppm	max 5 ppm
1,4-Dioxane	max 5 ppm	max 5 ppm
Sulfated ash	max 0.4 g/100 g	max 0.4 g/100 g
Volatile organic compounds	acc. to monographs	acc. to monographs

*The test methods are described in the current monographs USP-NF “Poloxamer”, JPE “Polyoxyethylene/Polyoxypropylene Glycol” and Ph. Eur. family monograph “Poloxamers”.*

## Solubility:

Both products are readily soluble in water and ethanol (95 %) and other mainly polar solvents. They are insoluble in ether, paraffin and fatty oils.

## Rheological properties:

Diluted aqueous solutions of both products show Newtonian flow properties that change to plastic flow properties. For **Lutrol®micro 68** it changes at a concentration of approx. 60% and for **Lutrol®micro 127** it changes at concentrations of above 15 %.

**Lutrol®micro F 127** is thermoreversible. **Lutrol®micro 68** shows a maximum viscosity at 60 to 75 °C whereas **Lutrol®micro 127** shows the maximum between 30 to 60 °C.




## The Application

**Lutrol®micro** grades are the excipients of choice for your dissolution problems with solid oral dosage forms.

**Lutrol®micro grades offer:**


- A high HLB value of > 20
- An excellent miscibility
- A wide acceptance of regulatory authorities



*Conventional nonionic surfactants are believed to act as a physical barrier to separate particles*

*Beyond CMC (critical micelle concentration), addition of more surfactant does not provide additional stabilization*

*Over time, aggregation can occur*




**Lutrol®micro grades can be used in several pharmaceutical application fields:**

- Dissolution enhancer for actives in tablets and capsules
- Polishing agent for film-coated tablets
- Lubricant for actives incompatible with magnesium stearate, e.g. ibuprofen
- Dispersing/wetting agent
- Water soluble lubricant e.g. for effervescent tablets

Due to the low average particle size that is in the same range as for active ingredients a direct blending with the active is possible with a low tendency of demixing.

**This allows a very cost- and time-saving production.**

*There is evidence that Lutrol®micro molecules do not micellize in the classical sense. Molecular aggregation is believed to occur over a broad concentration range, creating a “stacking effect” with a thicker shear plane.*





## The Formulation

*To show the good properties of Lutrol®micro grades a number of poorly soluble actives was tested.*


In this case a very easy manufacturing method was chosen to show only the influence of the solubilizer. Active ingredients and **Lutrol®micro** were blended for 10 minutes.

After mixing the material was filled into capsules and the dissolution profile was determined in an aqueous dissolution medium.


To give detectable values of the reference material (without solubilizer) the dissolution medium was mixed with ethanol.



**Lutrol® F**



**Lutrol® micro**






It can be shown that **L $\mu$ trol<sup>®</sup>micro 68** improves the dissolution of ibuprofen significantly. A very high release rate can be achieved by only blending both materials.



**Ibuprofen dissolution with and without L $\mu$ trol<sup>®</sup>micro**



- Ibuprofen 40 mg + 200 mg L $\mu$ trol<sup>®</sup>micro 68
- Ibuprofen 40 mg


*Ibuprofen crystals*

## The Formulation

The same behavior can be shown with the very poorly soluble active 17- $\beta$ -estradiol. The addition of **L $\mu$ trol<sup>®</sup>micro** increased the

drug release after 20 minutes from 15 % up to approx. 50 %. These results show the good solubilizing properties of **L $\mu$ trol<sup>®</sup>micro**.

### 17- $\beta$ -estradiol dissolution with and without **L $\mu$ trol<sup>®</sup>micro**



- 20 mg 17- $\beta$ -estradiol + 100 mg **L $\mu$ trol<sup>®</sup>micro 68**
- 20 mg 17- $\beta$ -estradiol






## The Summary

11

- Poloxamers can be used in controlled drug delivery systems leading to optimized bioavailability.
- The particle size of the micronized poloxamer is close to that of active ingredients. That perfects content uniformity.
- Micronized poloxamers can be used effectively in dry granulation and to minimize process related problems.

***Lptrol®micro grades are mostly applicable in solid dosage forms.***



Please write  
in block letters.  
Thank you.

Title/Name

Department

Company

Address

Town

Country

Telephone

Fax

e-Mail

Please  
apply  
postage  
stamp.

Reply card

**BASF**  
**Aktiengesellschaft**  
**G-MEP/ME – Li 554**  
**Attn: Dr. Hubertus Folttmann**  
**Carl-Bosch-Straße 64**  
**D-67117 Limburgerhof**  
**Germany**

## Fax Reply

Local contact:

or Headquarter Germany:

**+49-621-60-2 86 40**

Please complete, copy and fax to us,  
or detach the postcard and send it to us.

**Please send the  
following information.**

■ Technical information on  
**L $\mu$ trol<sup>®</sup> micro.**

■ Sample of  
**L $\mu$ trol<sup>®</sup> micro 68**, 0.5 kg.

■ Sample of  
**L $\mu$ trol<sup>®</sup> micro 127**, 0.5 kg.

■ Please contact me, I would  
like to know more about  
**L $\mu$ trol<sup>®</sup> micro.**

■ Technical information  
on Kollicoat<sup>®</sup> IR White.

■ Technical information  
on Kollicoat<sup>®</sup> Protect.

■ Technical information on  
Kollicoat<sup>®</sup> IR.

■ DVD "Pharmaceutical Excipients  
by BASF Fine Chemicals".

■ Newsletter "ExAct"  
(Excipients & Actives for Pharma).



**Please send the following information.**

- Technical information on **Lutrol®micro**.
- Sample of **Lutrol®micro 68**, 0.5 kg.
- Sample of **Lutrol®micro 127**, 0.5 kg.
- Please contact me, I would like to know more about **Lutrol®micro**.
- Technical information on Kollicoat® IR White.
- Technical information on Kollicoat® Protect.
- Technical information on Kollicoat® IR.
- DVD "Pharmaceutical Excipients by BASF Fine Chemicals".
- Newsletter "ExAct" (Excipients & Actives for Pharma).

**We look forward to answering any questions you may have. Please fill in the postcard, detach it and return it to the address overleaf.**

Please contact your local BASF company or one of the following regional centers:

#### **Asia**

**BASF Asia Pacific Regional HQ**  
Pharma Solutions  
Dr. Danilo Mercado  
BASF East Asia Regional Headquarters Ltd.  
45th Floor, Jardine House,  
No.1 Connaught Place,  
Central, Hong Kong  
Phone: +852-27311-588  
Fax.: +852-27349-638  
mercadd@  
basf-east-asia.com.hk

#### **Europe**

**BASF Aktiengesellschaft**  
MEE/HP – J 550  
Mr. Peter Hoffmann  
D-67056 Ludwigshafen  
Germany  
Phone: +49-621-60-7 69 28  
Fax: +49-621-60-7 69 40  
peter.01.hoffmann@  
basf.com

#### **North America**

**BASF Corporation**  
Pharma Solutions  
Mr. Javier Beeck  
100 Campus Drive  
Florham Park, NJ 07932  
USA  
Phone: +1-973-245-6393  
Fax: +1-973-245-6763  
javier.beeck@  
basf.com

#### **South America**

**BASF S.A.**  
Human Fine Chemicals  
Ms. Vanessa Occhipinti  
Estrada Samuel  
Aizemberg, 1707  
09851-550 São Bernardo  
do Campo – SP  
Brazil  
Phone: +55-11-4343-3311  
Fax: +55-11-4343-2255  
vanessa.occhipinti@  
basf.com

**[www.pharma-solutions.basf.com](http://www.pharma-solutions.basf.com)**

#### ■ Excipients

**Kollidon® grades**  
Group of povidone and copovidone products suitable mainly as tablet binders, crospovidone as tablet disintegrant and dissolution enhancer.

**Kollidon® SR**  
Matrix sustained release polymer.

**Ludipress® grades**  
Direct tableting aids for faster product development and speedier processing.

**Kollicoat® grades**  
Range of aqueous based film formers, cost efficient and ecological.

**Cremophor® grades and Solutol® HS 15**  
Range of different ethoxylated emulsifiers and solubilizers suitable for topical, oral and parenteral formulations.

**Soluphor® P**  
2-pyrrolidone.

**Lutrol® grades**  
Range of PEGs (Lutrol® E range) and poloxamers (Lutrol® F range) for a wide variety of pharmaceutical dosage forms.

#### ■ Actives

Ephedrines  
Pseudoephedrines  
Theophylline  
Caffeine  
Isotretinoin  
Tretinoin  
Ibuprofen  
Dobutamine  
Dopamine  
Isometheptene mucate  
Oxymetazoline  
PVP-Iodine  
Selegiline  
Xylometazoline

#### APIs by Orgamol

#### Vitamins

#### ■ Contract Manufacturing

#### ■ Value Added

## Pharma Solutions by BASF

BASF offers more than cGMP quality and supply safety: *technical expertise*. Our technical service is always at your side.

BASF Aktiengesellschaft  
Fine Chemicals Division  
Pharma Solutions  
67056 Ludwigshafen  
Germany

BASF wishes to create a prosperous and sustainable future with you as our customer – **and partner**.

■ Fax +49-621-60-2 86 40  
■ [pharma.solutions@basf.com](mailto:pharma.solutions@basf.com)  
■ [www.pharma-solutions.basf.com](http://www.pharma-solutions.basf.com)

 **BASF**

The Chemical Company

**BASF – the world's leading chemical company – can look back on well over 140 years of success and has attained an outstanding position as a reliable partner.**

**Our portfolio for the pharmaceutical industry comprises a comprehensive range of major and new active ingredients and excipient brands.**

